barybody

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A-Type Guide

A-Type | Trigger Point

Pointed shape for trigger points requiring precise, targeted massage

#Cervical Spine

#Trapezius

#Scalp

Cervical Spine (Neck)

STEP 1. Place the A-Type section facing upward and sit down.

STEP 2.

Lean back and lie down so that the two pressure points press against the center of the back of your neck.

STEP 3.

Slowly move your head left and right to release tension in the neck muscles.

Trapezius

STEP 1.

Position the A-Type so that the two pressure points sit on either side of the neck.

STEP 2.

Press firmly and massage from the neck down toward the trapezius muscle.

Scalp

Place the A-Type section gently on the scalp and apply steady pressure, massaging around the crown of the head.

B-Type Guide

B-Type | Large Muscle

Large, rounded shape designed to massage broad muscle groups with ease

#Thoracic Spine

#Lumbar

#Glutes

#Hip Joint

Thoracic (Upper Back) & Lumbar (Lower Back)

STEP 1. Place the B-Type section facing upward and sit down.

STEP 2. Lie back, rounded surface against the center of your back.

STEP 3. Slowly shift your body left and right to relax.

Lumbar (Lower Back)

STEP 1.

Lie face up and place the massager under your lower back.

STEP 2.

Position the lumbar area over the two pressure points and press to find the right spot.

STEP 3.

Slowly move your upper body left and right to massage.

STEP 4.

For deeper pressure, cross one leg over the other.

Scalp

Using the B-Type, press gently on the temples and massage in circular motions.

C-Type Guide

C-Type | Acupressure & Gua Sha

C-shaped pressure points designed for effective acupressure along the body's curves

#Lymph Nodes

#Calf

#Sole


Lymph Nodes

Position the C-Type below the collarbone and apply gentle circular pressure with moderate force.

Calf

Sit with knees bent and use the C-Type to press firmly along the calf, finding the right pressure points.

Sole

STEP 1.

Place the C-Type section facing upward on the floor.

STEP 2.

Rest one foot at a time on the pressure points and massage.

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